

COTTON



SLK-200-COTTON

MAGNETIC RESONANCE ANALYZER NMR



SIMULTANEOUS DETERMINATION

Oil
Moisture
Oleic & linoleic acid

· Fast and non-destructive measurement.

ACCURACY TABLE

MEASURING 27 cm ³ SEEDS				
	REPEATABILITY (Reliability 95%)	ACCURACY (Reliability 95%)	MEASUREMENT TIME	
OIL	± 0,06 %	± 1 %	4 Sec.	
MOISTURE	± 0,07 %	± 0,6 %	4 Sec.	
FATTY ACIDS (*)	± 0,5 %	± 2 %	20 Sec.	
PROTEIN	± 0,2 %	± 0,07 %	12 Sec.	

(*) OLEIC ACID, LINOLEIC ACID.

RN 20 km 24,4 | B° Causana · Lote 5 Mz.83
(CP 5101) Malagueño · Córdoba · Argentina
PH +54 3541 486 628/629 · info@spinlock.com.ar
www.spinlock.com.ar