

# OLIVE

## SLK-200-OLIVE

MAGNETIC RESONANCE ANALYZER NMR

**DRY AND HUMID SAMPLES**

· Olive, olive oil and pomace analysis  
**SIMULTANEOUS DETERMINATION**

**Oil**  
**Moisture**  
**Fatty acids**

· Fast and non-destructive measurement.



### ACCURACY TABLES

MEASURING 27 cm <sup>3</sup> of GROUND OLIVE (DRY)			
	REPEATABILITY (Reliability 95%)	ACCURACY (Reliability 95%)	MEASUREMENT TIME
OIL	± 0,06 %	± 1 %	4 Sec.
MOISTURE	± 0,07 %	± 0,6 %	4 Sec.
FATTY ACIDS (*)	± 0,50 %	± 2 %	20 Sec.

(\*) OLEIC ACID, LINOLEIC ACID.

MEASURING 27 cm <sup>3</sup> of GROUND OLIVE (WET)		
	ACCURACY (Reliability 95%)	MEASUREMENT TIME
OIL	± 1 %	2 Min.
MOISTURE	± 2 %	2 Min.

RN 20 km 24,4 | B° Causana · Lote 5 Mz.83  
(CP 5101) Malagueño · Córdoba · Argentina  
PH +54 3541 486 628/629 · info@spinlock.com.ar  
[www.spinlock.com.ar](http://www.spinlock.com.ar)